

## **SOP: Propagation of K562 (ATCC CCL-243)**

### **Information**

Name: K562

ATCC #: CCL-243

Tissue: bone marrow

Product Format: frozen

Morphology: lymphoblast

Culture Properties: suspension

Biosafety Level: 1

Disease: chronic myelogenous leukemia; 53-year-old female

ENCODE Number: HO2074

### **Materials List**

1. RPMI1640 (Life Technologies; Cat# 11875-150)
2. Heat Inactivated Fetal Bovine Serum (Life Technologies; Cat# 10082147)
3. Antibiotic-Antimycotic 100X (Life Technologies; Cat#15240-122)
4. T25, T75, T150, T182 culture flasks
5. Graduated pipets (1, 5, 10, 25, 50mL)
6. Freezing medium (Complete growth medium 95%; DMSO, 5%)
7. DMSO (Fisher; Cat#BP-231-100)
8. Cryovials (Sarstedt; Cat #72-694-006)
9. TC20 cell counter (Bio-Rad)
10. Counting Slides (Bio-Rad; Cat# 145-0011)
11. Microscope

### **Growth Medium for K562**

RPMI1640

10% FBS

Anti-Anti; 5 mL per 500 mL of culture medium

### **Procedure**

#### **A. Receipt of Frozen Cells and Starting Cell Culture**

- 1) Immediately place frozen cells in liquid nitrogen freezer storage until ready to culture.
- 2) When ready to start cell culture, quickly thaw ampoule in a 37°C water bath.
- 3) As soon as ice crystals disappear, swab outside surface of the ampoule with 70% ethanol, then transfer contents of ampoule into a 15 mL centrifuge tube with 9-10 mL of warm growth media.
- 4) Centrifuge at 1,500 rpm for 3 minutes.
- 5) Aspirate the medium. Suspend the cells in 10 mL complete growth medium and transfer the cells to a T25 cell culture flask.

6) Incubate at 5% CO<sub>2</sub>, 37 °C, humidified incubator. Keep the flask horizontally.

Note: These cells get into log phase in 5 to 7 days. Start counting the cells at day 3. When the total cell number reaches 0.7 to 0.8 X 10<sup>6</sup>/mL, split the culture to about 0.4 X 10<sup>6</sup> per mL fresh medium (you could use larger size flasks). From this point on the cells should double every 24 hours.

7) Record each subculture event as a passage.

## **B. Maintenance and Generation of Seed Stocks**

1) Following first or second passage after receipt of cells and with sufficient number of cells to continue maintenance and expansion. A small portion should be set aside as a seed stock. The cell pellet for the seed stock should be resuspended in freezing medium.

2) Cells in freezing medium are dispensed into cryovials (1-2 million cells per 1mL aliquot) and frozen in a -80°C cryo-freezing container overnight.

3) Cryovials are transferred the next day to liquid nitrogen freezer for long-term storage.

## **C. Harvest**

1) Passage cells until the desired number of cells is reached.

2) Remove cells from flasks as described above under “Sub-culture”.

3) Examine viability using Trypan blue staining.

Note: In our lab, we thaw a fresh vial of frozen cells, and grow up the cells for 4 to weeks.