

Supplementary Table S1. Primers used for reverse-transcriptase (RT) and real-time quantitative (RTQ) PCR

RT-PCR primers	Sense	Antisense
mGPR55	5'-GTGGTCTCCTTTCTCCCAGTG-3'	5'-TTAGCCCCTGGAGACCATGGTATCC-3'
mc-Fos	5'-AGCATGGGCTCTCCTGTCAA-3'	5'-GGGCTGCCAAAATAAACTCC-3'
GAPDH	5'-GGGAAGCTCACTGGCATGGCCTTCC-3'	5'-CATGTGGGCCATGAGGTCCACCAC-3'

RTQ-PCR primers	Sense	Antisense
hGPR55	5'-CTGCCTTGGTTCCACCATA-3'	5'-CCAGGATGCAGGTGAGTAAGA-3'
mGPR55	5'-TGGCCAGGCATCTTCAGT-3'	5'-CCAAGAGAAGTCCCCTTTCC-3'
mIL1 β	5'-TGTAATGAAAGACGGCACACC-3'	5'-TCTTCTTTGGGTATTGCTTGG-3'
mTNF α	5'-TGCCTATGTCTCAGCCTCTTC-3'	5'-GAGGCCATTTGGGAACCTTCT-3'
18S	5'-GCTCTAGAATTACCACAGTTATCCAA-3'	5'-AAATCAGTTATGGTTCCTTTGGTC-3'